



MDS3030K06+MD30-ADJSHLF2-XBAR

with Adjustable Overhead Shelf

M-SERIES TALL MEDICAL EQUIPMENT CART

with MD30-ADJSHLF2-XBAR Accessory

Tall, Standard Width, Six Drawer Monitor Cart with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

Accessories Include:

Adjustable Overhead Shelf (MD30-ADJSHLF2-XBAR)



Features:

- Sturdy, uni-body steel construction
- Choice of lock options: www.harloff.com/locking-options/
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Antimicrobial plastic top, removable
- Fully customizable drawer configurations
- Ergonomic side push handles can be mounted on either side
- 30" of vertical drawer space
- Cart is delivered assembled
- Latex free

M-series

Color shown is Navy. More than 20 additional colors are available.

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 43.75" H x 37.5" W x 22" D (111.1cm H x 95.3 cm W x 55.9 cm D) *without shelves*
- Cart weight: 174 lbs (79 kg)
- Approximate shipping weight: 213 lbs (96.6 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com