



**M-series**

**MDS3030K06+MD30-ANS**

## M-SERIES TALL ANESTHESIA CART

with MD30-ANS Anesthesia Package

### Tall, Standard Width, Six Drawer Anesthesia Cart with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

### MD30-ANS Accessories Include:



3-Tier Back Rail  
(MD30-RBR3)



(3) 12" Side Rails  
(SIDERAIL12)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



### Features:

- Sturdy, uni-body steel construction
- Choice of lock options: [www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Removable plastic top
- Fully customizable drawer configurations
- 30" of vertical drawer space
- Cart is delivered assembled
- Latex free

Color shown is Navy.  
More than 20 additional  
colors are available.

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**[medicalsales@harloff.com](mailto:medicalsales@harloff.com) • [www.harloff.com](http://www.harloff.com)**