

M-series MDS3030K06+MD30-ANS3

M-SERIES TALL ANESTHESIA CART

with MD30-ANS Anesthesia Package

Drawer Configuration:

121 cm W x 55.9 cm D) Cart weight: 224 lbs (101.6 kg)

Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm) 24.5" W x 17" D (62.2 x 43.2 cm) **Assembled Cart Dimensions:**

66.75" H x 47.625" W x 22" D (169.5 cm H x

Approximate shipping weight: 263 lbs (119.3 kg)

Approximate Shipping Dims: 56" H x 37" W x

26" D (142.2 cm x 94 cm x 66 cm) Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

MD30-ANS3 Accessories Include:




5-Compartment Tilt Bin w/Rail Clips



Plastic Waste Container (TILTBIN5RC) (MD-3GALWASTEDM) (MD-SHARPS-RC)



Sharps Container w/Wire Holder



Catheter Box w/ Rail Clip (CATHBOX18RC)



in the U.S.



3-Tier Back Rail (MD30-RBR3)



(3) 12" Side Rails (SIDERAIL12)



Large Utility Basket (BASKET3RC)



Divider with Trav



Color shown is Navy. More than 20 additional colors are



Tape and Label Dispenser w/Rail Clips (LABELDISP)



(1) Basic 6" Drawer Divider (MD30-DIV6-B)



(2) Standard 3"Drawer (MD30-TRAYDIV3-S)

Features:

- Sturdy, uni-body steel construction
- Choice of lock options: www.harloff.com/locking-options/
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling Tente® casters 1 steering caster and 2 braking casters
- Antimicrobial plastic top, removable
- Fully customizable drawer configurations
- 30" of vertical drawer space
- Cart is latex free and delivered assembled





View similar products online at

https://www.harloff.com/product-category/anesthesia-carts/

Toll Free: (800) 433-4064 • Ph: (719) 637-0300 medicalsales@harloff.com • www.harloff.com