



# MDS3030K16+MD30-AIRWAYPKG

## M-SERIES TALL TREATMENT CART with MD30-AIRWAYPKG Difficult Airway Package

### MD30-AIRWAYPKG Accessories Include:



Oxygen tank holder (O2HLDR)



(2) Utility hooks (UHOOKS2)



(1) 16" Side Rail (SIDERAIL16)



Adjustable Overhead Shelf (MD30-ADJSHLF2)



5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)



4-Scope Capacity Side Scope Cabinet (680350)

Color shown is Navy. More than 20 additional colors are available.

Engineered & Manufactured in the U.S.



### Features:

- Sturdy, uni-body steel construction
- Choice of lock options: [www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Aluminum side rail/push handle mounted on one side
- Industry best 5" quiet rolling casters - one steering caster and one braking casters
- Removable plastic top
- Fully customizable drawer configurations
- 30" of vertical drawer space
- Cart is latex free and delivered assembled

**M-series**

#### Drawer Configuration:

- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) & one 12" drawer (30.5 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
[medicalsales@harloff.com](mailto:medicalsales@harloff.com) • [www.harloff.com](http://www.harloff.com)