



# MDS3030KC4+MD30-DRW9-EH Equipment Holder

**M-SERIES STANDARD WIDTH CART**  
with 9" Equipment Holder, 3 Drawers and Key Lock

## M-Series Features:

- Sturdy, uni-body steel construction
- Choice of lock options:  
[www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Up to 30" of fully customizable vertical space is available for both drawer and equipment holder configurations
- 17" deep open equipment holding area can replace standard drawer areas in three possible height options
- 2" pass-through hole in the rear for cords and cables
- Full extension drawer slides for easy access
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Industry only fully curved protective bumper
- Durable scratch resistant hard baked powder coat finish
- Body has pre-threaded accessory mounting holes
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

## Colors:

Shown in Red, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).



# M-series

### Overall Cart Dimensions:

- Internal Drawer Dimensions: 24.5" W x 17" D (62.2 cm W x 43.2 cm D)
- Work Surface Height: 42.625" (108.3 cm)
- Cart dims: 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9cm D)
- Cart net weight: 131 lbs (59.4 kg)
- Approx shipping weight: 168 lbs (76.2 kg)
- Approx shipping dims: 56" H x 39" W x 26" D (142.2 cm H x 99 cm W x 66 cm D)

### Drawer Configuration:

- One MD30-DRW9-EH equipment holder with an 8" H x 24.625" W x 17" D (20.3 cm H x 62.5 cm W x 43.2cm D) opening, one 3" drawer (7.6 cm) and two 9" drawers (22.9 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Designed and built to meet the needs of today's healthcare market. Replace 9", 12 or 18" of vertical drawer space at the top or bottom of your M-Series cart with a doored equipment holder, perfect for storing electronics or large medical equipment. M-Series carts offer you choices.

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300