

Malignant Hyperthermia Cart with Medical Grade Refrigerator/ Model: MH5300B

Three Drawer with Breakaway Lock



Eighteen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

- * NOTE: Refrigerator units are subject to change without notice.
- * NOTE: Must be used with 120 Volt 50/60 Hertz electrical power and NEMA 5-15 Wall Socket.

Construction:

- 20-gauge powder-coated steel with built-in Follett® REF2 refrigerator

Features:

- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where they are most needed
- Lab-quality Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawers and refrigerator, (refrigerator can be keylocked)
- Cart ships ready to use with wheels installed
- Hospital grade plug and cord included

Accessories Include:

- Waste container with mounting bracket
- Large wire basket
- Two large chart/document holders

Optional Accessories Include:

- Side storage cabinets with key locking door and two adjustable shelves for left side (MH-SIDECABL) and/or right side (MH-SIDECABR) for additional supply storage (each side cabinet replaces existing side mounted accessories).

Drawer Configuration:

- (3) – 3" (7.6 cm)
- (1) – Compartment with pre-installed refrigerator

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall Cart Dimensions:

- Cart dims: 51.5 H x 37.5" W x 26"D (130.8 cm H x 95.3 cm W x 66 cm D)
- Cart net weight: 188 pounds (85.3 kg)
- Approximate shipping weight: 227 pounds (103.2 kg)
- Approximate shipping dims: 62" H x 44" W x 30" D (157.5cm H x 111.8cm W x 76.2cm D)