



MALIGNANT HYPERTHERMIA

MH5300K

MALIGNANT HYPERTHERMIA CART with 1.8 Cubic Feet Medical Grade Refrigerator

Features:

- 20-gauge powder-coated steel construction w/built-in Follett® REF2 refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- Highest quality in the industry Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Key locks for drawer and refrigerator
- Cart ships ready to use with wheels installed
- Antimicrobial plastic top, removable

Colors:

Shown in Hammertone Blue, but available in over twenty color options. View them online: www.harloff.com/finish-options/



Malignant Hyperthermia Cart, Three Drawer with Key Lock

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Overall Cart Dimensions:

- Internal Drawer Dimensions: 23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 51.5" H x 37.5" W x 26" D (130.8 cm H x 95.3cm W x 66cm D)
- Cart net weight: 188 lbs (85.3 kg)
- Approximate shipping weight: 259 lbs (117.5 kg)
- Approximate shipping dims: 59" H x 33" W x 28" D (150 cm H x 83.8 cm W x 71.2 cm D)

Drawer Configuration:

- (3) 3" (7.6 cm) and (1) Compartment with pre-installed refrigerator

Accessories Included:



Large wire basket



Waste container & mounting bracket



Two large chart/document holders



Engineered & Manufactured in the U.S.

OPTIONAL:

Side storage cabinets with key locking door and two adjustable shelves for additional supply storage (each side cabinet replaces existing side mounted accessories).



View similar products online at <https://www.harloff.com/malignant-hyperthermia-cart/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300