

## Malignant Hyperthermia Cart with Medical Grade Refrigerator/ Model: MH5300K

Three Drawer with Key Lock



### Eighteen Available Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

- \* NOTE: Refrigerator units are subject to change without notice.
- \* NOTE: Must be used with 120 Volt 50/60 Hertz electrical power and NEMA 5-15 Wall Socket.

### Construction:

- 20-gauge powder-coated steel with built-in Follett® REF2 refrigerator

### Features:

- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where they are most needed
- Lab-quality Follett® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Key locks for drawers and refrigerator
- Cart ships ready to use with wheels installed
- Hospital grade plug and cord included

### Accessories Include:

- Waste container with mounting bracket
- Large wire basket
- Two large chart/document holders

### Optional Accessories Include:

- Side storage cabinets with key locking door and two adjustable shelves for left side (MH-SIDECABL) and/or right side (MH-SIDECABR) for additional supply storage (each side cabinet replaces existing side mounted accessories).

### Drawer Configuration:

- (3) – 3" (7.6 cm)
- (1) – Compartment with pre-installed refrigerator

### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

### Overall Cart Dimensions:

- Cart dims: 51.5 H x 37.5" W x 26"D (130.8 cm H x 95.3 cm W x 66 cm D)
- Cart net weight: 188 pounds (85.3 kg)
- Approximate shipping weight: 227 pounds (103.2 kg)
- Approximate shipping dims: 62" H x 44" W x 30" D (157.5cm H x 111.8cm W x 76.2cm D)