



*Wood Style*

**MPA3024E04-AN**

**ALUMINUM**

## **A-SERIES SHORT WOOD STYLE CART**

with Five Drawers and Electronic Pushbutton Lock

### **Why Wood Style Drawer Fronts?**

*Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.*

### **Features:**

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Basic electronic pushbutton lock with key lock override and factory set user/supervisor codes
- Cabinet may be programmed in the field
- Heat-sealed durable scratch resistant powder coat finish on body, inside and out
- Pre-threaded mounting holes for accessories found at [harloff.com/product-category/accessories/](http://harloff.com/product-category/accessories/)
- Industry only fully curved protective bumper guards your cart and walls from damage
- Removable plastic top, easily replaceable
- Ergonomic side push handles can be mounted on the left or right
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- 24" of vertical drawer space with fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters
  - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled



**Body shown in Beige. More than 20 additional colors are available. Drawers shown in Wood Vinyl American Natural.**

Three standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.



**Engineered & Manufactured in the U.S.**

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of A-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff aluminum cart.

#### **Drawer Configuration:**

- One 3" drawer (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



**View similar products online at**

<https://www.harloff.com/product-category/a-series-medical-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
www.harloff.com • Ph: (719) 637-0300**