



MPA3030B09PED1

ALUMINUM

A-SERIES TALL PEDIATRIC CART with Nine Drawers and a Breakaway Lock

Designed for treating pediatric emergencies by using a multi-colored drawer labeling system comparable to the Broselow color coding method

Features:

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on the left or right
- Heat sealed powder coat steel full length breakaway lock bar with 100 breakaway seals
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space



Color shown with a Light Gray body and color coded drawers. More than 20 additional body colors are available.



Tall, Standard Width, Nine Drawer Pediatric Cart with a Breakaway Drawer Lock

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.

Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer

Engineered & Manufactured in the U.S.

Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Writing surface: 42.625" (108.2 cm)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



View similar products online at <https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300