



Wood Style

MPA3030BC6-CM

ALUMINUM

## A-SERIES TALL WOOD STYLE CART

with Six Drawers and a Breakaway Lock

### Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.



Engineered & Manufactured in the U.S.

### Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Heat sealed powder coat steel breakaway lock bar with a pack of 100 break-away seals
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish on body, inside and out
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, removable and easily replaceable
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Ergonomic side push handles can be mounted on the left or right
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- 30" of vertical drawer space with fully customizable drawer configurations at no extra charge
- Pre-threaded mounting holes for accessories found at [harloff.com/product-category/accessories/](http://harloff.com/product-category/accessories/)
- Cart is delivered assembled

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of A-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff aluminum cart.



Body shown in Beige. More than 20 additional colors are available. Drawers shown in Wood Vinyl Cherry Mahogany.

Two standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.

#### Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300