



MPA3030E06+MD30-ANS3

ALUMINUM

A-SERIES TALL ANESTHESIA CART

w/Six Drawers, E-Lock & MD30-ANS3 Anesthesia Pkg

MD30-ANS3

Accessories Include:



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



Plastic Waste
Container
(MD-3GALWASTEDM)



Sharps Container
w/Wire Holder
(MD-SHARPS-RC)



Catheter Box w/
Rail Clip
(CATHBOX18RC)

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



3-Tier Back Rail
(MD30-RBR3)



(3) 12" Side Rails
(SIDERAIL12)



Large Utility Basket
(BASKET3RC)



Engineered &
Manufactured
in the U.S.



Tape and Label
Dispenser w/Rail
Clips (LABELDISP)



(1) Basic 6"
Drawer Divider
(MD30-DIV6-B)



(2) Standard 3"Drawer
Divider with Tray
(MD30-TRAYDIV3-S)

Color shown is
Navy. More than
20 additional
colors are
available.

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction
- Heat-sealed durable scratch resistant powder coat finish
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumpers
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top, easily replaceable
- Fully customizable drawer configurations at no extra charge
- Ergonomic side push handles, mounted on the left or right
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- 30" of vertical drawer space



HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/acute-care-carts/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com