

**MPA3030E16** with  
MD30-SSFLTOP accessory


## DELIVERY CART

with Stainless Steel Top + Two Drop Shelves

*The cart top is designed to create a large workspace ideal for a baby scales, pediatric supplies and other equipment*

### Features:

- Curved drawer pulls, durable and elegant
- Stainless steel top with stainless steel drop shelves on both sides
- The large top workspace (58.5" in length when shelves are fully extended) is ideal for holding a baby scale and other large equipment
- Basic pushbutton e-lock with key lock override
- Cabinet may be programmed in the field
- Factory set user and supervisor code
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters - 1 steering caster and 2 braking casters
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Latex free cart is delivered fully assembled
- 30" of vertical drawer space
- Drawer Configuration: Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown with a Light Blue body and drawers. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

**Tall, Standard Width Delivery Cart with Stainless Steel Top and Side Flyout Shelves, Keyless Pushbutton E-Lock**

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



View similar products online at

<https://www.harloff.com/product-category/a-series-medical-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**