



Wood Style

MPA3030EC6-CM

ALUMINUM

A-SERIES TALL WOOD STYLE CART with Six Drawers and Electronic Pushbutton Lock

Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.



Engineered & Manufactured in the U.S.

Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Fully welded uni-body aluminum construction, highest quality in the industry
- Basic electronic pushbutton lock with key lock override and factory set user/supervisor codes
- Cabinet may be programmed in the field
- Curved drawer pulls, durable and elegant
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, removable and easily replaceable
- Pre-threaded mounting holes for accessories found at harloff.com/product-category/accessories/
- Ergonomic side push handles can be mounted on the left or right
- Fully customizable drawer configurations at no extra charge
- Cart is delivered assembled
- 30" of vertical drawer space
- Industry only fully curved protective bumper guards your cart and walls from damage
- Heat-sealed durable scratch resistant powder coat finish on body, inside and out
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters



Body shown in White. More than 20 additional colors are available. Drawers shown in Wood Vinyl Cherry Mahogany.



Two standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of A-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff aluminum cart.

Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)

Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



View similar products online at <https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300