



MPA3030K06+MD30-ANS

ALUMINUM

A-SERIES STANDARD WIDTH CART w/Six Drawers, Key Lock & MD30-ANS Anesthesia Pkg

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.

MD30-ANS Accessories Include:



3-Tier Back Rail
(MD30-RBR3)



(3) 12" Side Rails
(SIDERAIL12)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



**Engineered &
Manufactured
in the U.S.**

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top, easily replaceable
- Fully customizable drawer configurations at no extra charge
- Ergonomic side push handles, mounted on the left or right
- Industry best 5" quiet rolling casters
 - 1 steeringcaster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space

**Color shown is Navy.
More than 20 additional
colors are available.**



Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com