



MPA3030K06+MDPW

**ALUMINUM**

## A-SERIES TALL COMFORT CART with Six Drawers and a Standard Key Lock

*Comfort carts provide tangible and psychological benefits to patients by distributing items for fun, relaxation and mental health rather than clinical needs*

### Features:

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction with heat sealed, durable powder coat finish
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Fully curved protective bumper
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, removable/replaceable
- Ergonomic side push handles
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Unmatched selection of alternative lock options
- 30" of vertical drawer space
- Cart is delivered assembled

*Cart options include full wrap or half wrap (shown to the right) with the exterior body covered and your choice of available drawer colors.*

**Wraps available for all steel M-Series & aluminum A-Series carts!**

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function.



**Cart is fully wrapped with an Alphabet theme. Color shown is Light Blue. More than 20 additional body/drawer colors are available.**



**Engineered & Manufactured in the U.S.**

#### Drawer Configuration:

- Three 3" (7.6 cm) drawers, two 6" (15.2 cm) drawers and one 9" (22.8 cm) drawer
- 24.5" W x 17" D (62.2 x 43.2 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (1111 cm H x 93.3 cm W x 55.9 cm D)
- Writing surface: 42.625" (108.2 cm)
- Approximate shipping dims: 56" H x 39" W x 26" D (142.2 cm x 99 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Top Image by Freepik



View similar products online at <https://www.harloff.com/product-category/pediatric-resuscitation-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**



Hawaiian Flora



Sunflowers



Evergreen



Balloons



Pawprints



Camouflage



Ambulance



Space



Snowflakes



Artistic Mountains



Angular Peaks



Alphabet



Dinosaurs



Bug Life



Handprints



Custom



Sports



Numbers

*Want something different than one of our pre-selected designs? Supply your own image!*

Our cart wraps can use repeatable patterns, like those shown above, or single images like hospital logos, photographs or virtually any other image with a high resolution\*. Call Harloff to speak to a representative about your requirements to customize your cart to your specifications.

\*Scalable vector files (.ai or .eps) work best, pixel based images (.psd, .tif, etc) may be acceptable if they have a high enough digital resolution (300 dpi minimum).