



MPA3030K16+MD30-AIRWAYPKG

ALUMINUM

A-SERIES TALL ANESTHESIA CART

6 Drawer Key Lock, MD30-AIRWAYPKG Anesthesia Pkg


MD30-AIRWAYPKG Accessories Include:



Oxygen tank
holder (O2HLDR)



(2) Utility hooks
(UHOOKS2)



(1) 16" Side Rail
(SIDERAIL16)



Adjustable
Overhead Shelf
(MD30-ADJSHLF2)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



4-Scope Capacity
Side Scope
Cabinet (680350)

Drawer Configuration:

- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Color shown is
Navy. More than 20
additional colors
are available.

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Ergonomic side push handles, mounted on the left or right
- Removable plastic top, easily replaceable
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters
 - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space



View similar products online at
<https://www.harloff.com/acute-care-carts/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com



Engineered &
Manufactured
in the U.S.