

MPA3030K16+MD30-AIRWAYPKG

aluminum

A-SERIES TALL ANESTHESIA CART

6 Drawer Key Lock, MD30-AIRWAYPKG Anesthesia Pkg

MD30-AIRWAYPKG Accessories Include:





Oxygen tank holder (O2HLDR)



Adjustable

(2) Utility hooks (UHOOKS2)



5-Compartment Overhead Shelf Tilt Bin w/Rail Clips (MD30-ADJSHLF2) (TILTBIN5RC)

- · Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- · Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to prethreaded holes with a screwdriver
- · Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Ergonomic side push handles, mounted on the left or right
- Antimicrobial plastic top, removable and easily replaceable
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space



Color shown is Navy. More than 20 additional colors are available.

Toll Free: (800) 433-4064 • Ph: (719) 637-0300

medicalsales@harloff.com • www.harloff.com





Engineered & Manufactured in the U.S.



Side Scope Cabinet (680350)

(1) 16" Side Rail (SIDERAIL16)

66.75" H x 37.5" W x 22" D (169.5 cm H x 95.3 cm W x 55.9 cm D) Cart weight: 198 lbs (89.8 kg) without accessories

Drawer Configuration:

Approximate shipping weight: 237 lbs (107.5 kg)

24.5" W x 17" D (62.2 x 43.2 cm) Assembled Cart Dimensions:

Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts