



MRI-SAFE

MR7B

MR-CONDITIONAL EMERGENCY CART with Seven Drawers and a Breakaway Lock

Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- Breakaway locking system (with optional multiple lock bars)
- Pack of 100 breakaway seals
- Protective corner bumpers
- Pull-out shelf
- Antimicrobial plastic top with integrated push handles, removable
- Latex free cart is delivered fully assembled with pre-threaded accessory mounting holes
- Durable, heat sealed powder coat finish



Engineered & Manufactured in the U.S.



Color shown is White with Red drawers. More than 20 additional colors are available.

Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)
- One 4" drawer (10.2 cm), four 3.25" drawers (8.3 cm), one 6.5" drawer (16.5 cm) and one 9.75" drawer (24.8 cm)

Carts Specially Made for Magnetic Resonance Imaging Suites

Harloff offers a line of carts that meet American Society for Testing and Materials (ASTM) International Designation of MR-Conditional for use in a 3 Tesla or less environment of the Magnetic Resonance Imaging Suite (testing date 1-13-2010). Whether you need it for Anesthesia or Emergency, Harloff has the right MR-Conditional cart for your needs.

Assembled Cart Dimensions:

- Writing surface: 43.5" (110.5cm)
- 46" H x 31.8" W x 23.75" D (116.8 cm H x 80.8 cm W x 60.3 cm D)
- Cart net weight: 110 lbs (49.9 kg)
- Approximate shipping weight: 151 lbs (68.5 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7cm H x 94cm W x 66cm D)

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at <https://www.harloff.com/product-category/mri-safe-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300