

MR-Conditional Line/ Model: MR7B

Aluminum Seven Drawer Anesthesia Cart, Breakaway Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Non-magnetic or "weakly magnetic" materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- Painted aluminum breakaway locking system
- Breakaway lock seals, numbered, yellow, pkg-100
- Protective corner bumpers
- Painted steel pull-out shelf
- Antimicrobial plastic top with integrated push handles, removable
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (4) - 3.25" (8.3 cm)
- (2) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Internal Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)

Overall Cart Dimensions:

- 45.75" H x 30.5" W x 23.75" D (116.2 cm H x 77.5 cm W x 60.3 cm D)
- Writing surface: 43.5" (110.5 cm)
- Cart net weight: 110 pounds (49.9 kg)
- Approximate shipping weight: 151 pounds (68.5 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7 cm H x 94 cm W x 66 cm D)