



# MRI-SAFE

## MRN6B

### MR-CONDITIONAL EMERGENCY CART with Six Drawers and a Breakaway Lock

#### Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- Breakaway locking system (with optional multiple lock bars)
- Pack of 100 breakaway seals
- Protective corner bumpers
- Antimicrobial plastic top with integrated push handles, removable
- Latex free cart is delivered fully assembled with pre-threaded accessory mounting holes
- Durable, heat sealed powder coat finish


#### Drawer Configuration:

- 13.75" W x 17" D (34.9 cm x 43.2 cm)
- One 4" drawer (10.2 cm), two 3.25" drawers (8.3 cm), two 6.5" drawers (16.5 cm) and one 9.75" drawer (24.8 cm)



Engineered & Manufactured in the U.S.

Color shown with a White body and Red drawers. More than 20 additional colors are available.



#### Carts Specially Made for Magnetic Resonance Imaging Suites

Harloff offers a line of carts that meet American Society for Testing and Materials (ASTM) International Designation of MR-Conditional for use in a 3 Tesla or less environment of the Magnetic Resonance Imaging Suite (testing date 1-13-2010). Whether you need it for Anesthesia or Emergency, Harloff has the right MR-Conditional cart for your needs.

#### Assembled Cart Dimensions:

- Writing surface: 44" (111.7 cm)
- 46" H x 21.5" W x 23.75" D (116.8 cm H x 54.6 cm W x 60.3 cm D)
- Cart net weight: 85 lbs (38.6 kg)
- Approximate shipping weight: 123 lbs (55.79 kg)
- Approximate shipping dims: 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)

**Warranty:** 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/mri-safe-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300