



- 1680 Box Capacity
- Mobile Body
- Glass Doors
- Standard Key Lock

**MEDSTOR<sup>MAX</sup>**

**MS-IOL1680**

## INTRAOCULAR LENS CART

### Features:

- Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications
- Integrated, recessed side push handles so cabinet can stack flush with walls or other cabinets — we don't use easily misplaced removable handles or bulky exterior hardware.
- Ball-bearing full extension drawer slides
- Features a true 24" interior cabinet depth for storage of up to 1680 IOL boxes
- Key lock with two keys for lockable, secure storage of supplies
- Eight pull-out drawer shelves with adjustable dividers — ideal for medical supplies used in surgery centers and operating rooms
- Cabinet is delivered fully assembled
- Double doors can swing up to 270 degrees for easy access
- Glass windows allow for visibility of supplies and inventory
- Industry best 5" quiet rolling Tente® casters for superior mobility - one steering and one braking caster



**Engineered & Manufactured in the U.S.**



#### Assembled Cart Dimensions:

- Cart net weight: 500 pounds (226.8 kg)
- 76" H x 44" W x 27.8" D (193cm H x 111.8cm W x 70.6cm D)
- Approximate shipping weight: 590 lbs (267.6 kg)
- Approximate Shipping Dims: 88" H x 48" W x 39" D (223.5 cm x 122 cm x 99.1 cm)

**Warranty:** 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Optional  
dimmable LED lighting  
is also available!**

### Mobile IntraOcular Lens Storage Cart, Full Height Cabinet, Eight Drawers, Glass Doors, Standard Key Lock

Attractive, secure storage for IntraOcular lenses. Provides mobility and is compatible with multiple brands and sizes of IOL boxes including Alcon, Bausch and Lomb, Rayner, and Abbot Medical Optics.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/intraocular-lens-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**