

Suture Storage and Transport Cart, Key Locking

#MS-SUTURE2-K

Designed for facilities requiring secure but visible storage and transport of sutures.



Tray not included as standard

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Color shown is sand

Assembled Cart Dimensions:

- 52.5" H x 33.5" W x 25"D
- Cart net weight: 243 pounds

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- 20-gauge powder-coated steel, STEEL FRAMED ACRYLIC DOORS

Features:

- Double sided with bottom pass-through
- Clear acrylic doors are steel-framed for durability with pull handles
- Clear doors allow inspection of contents
- 5 Adjustable shelves on each side
- Key lock for doors
- 5" Casters, all full swivel, 1 directional/tracking, one braking
- Cart ships assembled with wheels installed
- Durable, hard-baked powder coat finish available in more than twenty colors

Optional Trays:

- Bottom pass through can hold Harloff MedStorMax trays
 - 1 – 8" (81032S) or
 - 1 – 4" (81031S) and 1-2" or
 - 3 – 2" (81030S)

Choose Your Colors:

<https://www.harloff.com/finish-options/>

Choose Your Baskets/Trays:

<https://www.harloff.com/medical-storage-accessories/>