



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Double Column Short Mobile Medical Storage Cabinet with Key Locking Tempered Glass Double Doors #MS4140DR

The MS4140DR is a short, double column medical storage cart designed for quick access to medical instruments and supplies, featuring the most innovative modular tray system on the market.

The MS4140DR can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Standard Features:

- Plastic side panels with 19 insertion opportunities per column, to hold modular exchange trays and baskets,
- Key locking tempered glass double doors
- Painted steel top with integrated push handles
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand
- Twelve year warranty against manufactures defect

Advantages and Benefits:

- Maximize space and minimize costs using 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage, reducing capital purchase cost.
- Pull-out-stoppers available for trays, allowing trays to be viewed at an angle without falling out of the cart.
- Catheter shelf can be added

Assembled Cart Dimensions:

- 42.6" H x 40.1" W x 29.25" D
(108.2 cm H x 101.9 cm W x 74.3 cm D)
- Cart net weight: 193 pounds (87.5 kg)
- Approximate shipping weight: 330 lbs (149.7 kg)