

Preconfigured Double Column Medical Storage Cart with Key Locking Clear Panel Double Doors, #MS8140DR-A

The MS8140DR-A has key locking clear panel double doors, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8140DR-A can be used in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



Assembled Cart Dimensions:

- 76" H x 42" W x 29.75" D
(193.04 cm H x 106.68 cm W x 75.57 cm D)
- Cart net weight: 396 pounds (180 kg)
- Approximate shipping weight: 495 pounds (225 kg)
- Approximate shipping dims: 86" H x 48" W x 39" D
(218.4 cm H x 122 cm W x 99.1 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Offers 35% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 5" and 8" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique "clip-on" locking design to withstand high loads.
- Unique design of plastic side panels allows modular trays to be extended and viewed at an angle without falling out of the cart. (pull-out stoppers required)
- Catheter slide shelf can be added featuring three cath slides, and eight hooks per slide.

Standard Features:

- Two inch modular trays: Total of ten trays that are 23.50 x 15.75 x 2 inches
- Four inch modular trays: Total of eight trays that are 23.50 x 15.75 x 4 inches
- Eight inch modular trays: Total of six trays that are 23.50 x 15.75 x 8 inches
- Plastic side panels form a modular system consisting of 94 insertion opportunities (47 per column) to hold modular exchange trays and baskets.
- Key locking clear glass panel double doors
- Integrated push handles are flush with side of cart
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Ships fully assembled

Colors:

More than twenty color options are available.

Choose Your Colors Online:

<https://www.harloff.com/finish-options/>