



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Preconfigured Double Column Catheter Procedure Cart with Clear Panel Doors #MS8140DRCATH-3

The MS8140DRCATH-3 is a double column Catheter Procedure Cart with tempered glass double doors, specifically designed for use in a cath lab to support the catheterization procedure.

The MS8140DRCATH-3 comes preconfigured as follows: left column, top to bottom 1 catheter slide shelf with 3 slides; right column, top to bottom, 3ea 1/2" trays, 5ea 4" wired baskets, 4ea 7" wired baskets



Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- The significant amount of catheter sizes and packaging techniques has made organizing and managing this critical area an extremely difficult task. Harloff offers a solution with a high density Catheter Storage System, which provides quick access and retrieval of catheters and related supplies.
- Equipped with two adjustable Cath Slide Shelves, consisting of three full extension catheter slides, each slide has eight hooks, and each hook can hold up to 12 catheter packages.
- Each Cath Slide features a pair of ball bearing slides to withstand the pull force of removing catheters from hooks; full extension brings the catheters into full view for quick access and visual inventory of remaining stock.
- Cath Slide Shelf can be mounted in lower position for storage above.
- The MS8140DRCATH-3 stores long stent boxes, balloons, guidewires, and other supplies efficiently for easy access with our specially-engineered catheter and interior side panels. The maximum package length is 61 inches.

Standard Features:

- Catheter Slide Shelf, featuring three slides and eight hooks per slide, can be mounted in a lower position to create storage above for over-sized items.
- Cath slides consist of a pair of ball bearing drawer slides to withstand the pull force of removing catheters from hooks.
- Key locking tempered glass double doors
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty

Assembled Cart Dimensions:

- 76" H x 42" W x 29.6" D
(193 cm H x 106.7 cm W x 75.2 cm D)
- Cart net weight: 390 pounds (176.9 kg)
- Approximate shipping weight: 520 pounds (235.9 kg)