



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Triple Column Mobile Medical Storage Cabinet with Key Locking Roll-up Tambour Door MS8160CATH-3

The MS8160CATH-3 is a double column medical storage cabinet secured by a key locking roll-up tambour door. Utilizing the most innovative modular tray system on the market, this mobile cabinet can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.

Preconfigured as follows: left column, 3 ea MS-SH shelves, 2 ea 2" trays, 2 ea 4" trays, 2 ea 5" wired baskets, 1 ea 8" wired baskets; middle column, 3 ea MS-SH shelves, 4 ea 2" trays, 3 ea 4" trays, 2 ea 8" trays; right column, 1 ea CATH-3 catheter slide shelf



Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Maximize space and minimize costs using 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage, reducing capital purchase cost.
- Pull-out-stoppers available for trays, allowing trays to be viewed at an angle without falling out of the cart.

Standard Features:

- Plastic side panels with 47 insertion opportunities per column, to hold modular exchange trays and baskets,
- Full extension catheter holders with five pairs of hooks
- Key locking roll up tambour door
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand
- Twelve year warranty against manufactures defect

Assembled Cart Dimensions:

- 80.5" H x 59.6" W x 27.9" D
(204.5 cm H x 151.4 cm W x 70.8 cm D)
- Cart net weight: 400 pounds (181 kg)
- Approximate shipping weight: 590 pounds (267.2 kg)