



650 Ford Street  
Colorado Springs, CO 80915  
Toll: (800)433-4064 \* Ph: (719)637-0300  
Fax: (719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Triple Column Medical Storage Cart with Tempered Glass Doors #MS8160DR

Triple column cart with tempered glass doors, and plastic interior side panels of 141 slide tracks (47 per column) to hold various exchange trays and wired baskets\*

The MS8160DR can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be extended and viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains eight hooks.

### Standard Features:

- Plastic side panels form a modular system consisting of 141 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Catheter Shelf Holders can be added to any column.
- Integrated push handles are flush with cart side, positioned vertically for better control
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty

### Assembled cart dimensions:

- 76" H x 62.25" W x 27.9" D  
(193 cm H x 158.1 cm W x 70.8 cm D)
- Cart net weight: 474 pounds ( 215 kg)
- Approximate shipping weight: 599 pounds (271.7 kg)

\* Not Included