



- Powder Coat Steel
- Mobile Body
- Half Height
- Glass Door
- Standard Key Lock

# MEDSTOR<sup>MAX</sup>

## MSPM42-00GK


### DOUBLE COLUMN STORAGE CABINET

#### Features:

- Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications
- Largest selection of locking options in the industry including badge entry and WiFi-capable locks
- Improved, better refined aesthetics featuring sleek door handles and hinged doors
- Stacks flush with walls or other cabinets — we don't use easily misplaced removable handles or bulky exterior hardware
- Doors can swing up to 270 degrees, great for leaving cabinets open during busy days.
- Features a true 24" interior cabinet depth able to store boxes two feet deep
- Customize your cabinet with any combo of trays, baskets, shelves or accessories/packages
- Interior side panels allow 19 insertion points within each column creating a system of multiple organization possibilities
- Modular trays can be pulled out at an angle, and "clip-on" locking safety features help minimize the risk of falling trays
- Optional dimming LED lights toggle on/off when the door is open, ideal for use in surgical areas for better procedural vision
- Industry best 5" quiet rolling Tente® casters for superior mobility - 1 steering and 1 braking caster
- Cabinet is delivered fully assembled



Engineered & Manufactured in the U.S.



#### Assembled Cart Dimensions:

- 40.75" H x 40.1" W x 29.25" D  
(103.5 cm H x 101.9 cm W x 74.3 cm D)
- Cart net weight: 193 pounds (87.5 kg)
- Approximate shipping weight: 336 lbs (152.4 kg)
- Approximate Shipping Dims: 62" H x 49" W x 39" D  
(157.5 cm x 124.5 cm x 99.1 cm)

**Warranty:** 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Shown with  
(6) 81030S 2"  
trays, (4) 81031S  
4" trays (2)  
81032S 8" trays,  
sold separately

Shown in Yellow. Also available in more than 20 color options or Stainless Steel.

#### Mobile Cabinet, Half Height, Double Column, Full Panels, Glass Door, Standard Key Lock

For mobile or fixed storage, our versatile medical storage solutions provide high density and secure storage for valuable medical supplies used in surgical procedures. These med carts are designed around the most innovative modular tray system on the market and feature your choice of door style for security and ease of access.



View similar products online at

<https://www.harloff.com/product-category/medstor-max-medical-storage/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300