



- Powder Coat Steel
- Mobile Body
- Full Height
- Glass Door
- Standard Key Lock

MEDSTOR^{MAX}

MSPM84-00GK

QUAD COLUMN STORAGE CABINET

Features:

- Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications
- Largest selection of locking options in the industry including badge entry and WiFi-capable locks
- Improved, better refined aesthetics featuring sleek door handles and hinged doors
- Stacks flush with walls or other cabinets — we don't use easily misplaced removable handles or bulky exterior hardware
- Doors can swing up to 270 degrees, great for leaving cabinets open during busy days.
- Features a true 24" interior cabinet depth able to store boxes two feet deep
- Customize your cabinet with any combo of trays, baskets, shelves or accessories/packages
- Interior side panels allow 49 insertion points within each column creating a system of multiple organization possibilities
- Modular trays can be pulled out at an angle, and "clip-on" locking safety features help minimize the risk of falling trays
- Optional dimming LED lights toggle on/off when the door is open, ideal for use in surgical areas for better procedural vision
- Industry best 5" quiet rolling Tente® casters for superior mobility - 1 steering and 1 braking caster
- Cabinet is delivered fully assembled



Shown in Light Gray. Also available in more than 20 color options or Stainless Steel.

Shown with trays and wire baskets, sold separately

Assembled Cart Dimensions:

- 77.75" H x 82.125" W x 30.25" D (192.4cm H x 208.6cm W x 76.8cm D)
- Cart net weight: 980 lbs (444.5 kg)
- Approximate shipping weight: 1000 lbs (453.6 kg)
- Approximate shipping dims: 88" H x 96" W x 39" D (223.5 cm H x 243.8 cm W x 99.1 cm D)

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

For mobile or fixed storage, our versatile medical storage solutions provide high density and secure storage for valuable medical supplies used in surgical procedures.



View similar products online at
<https://www.harloff.com/product-category/medstor-max-medical-storage/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300



Engineered & Manufactured in the U.S.