



MWA2427K01-21-DG with
MD24-SSFLTOP accessory

DELIVERY CART

with Stainless Steel Top + Two Drop Shelves

The cart top is designed to create a large workspace ideal for a baby scales, pediatric supplies and other equipment

Features:

- **The wood style cart wrap gives the cart a furniture feel and is available in multiple finish options**
- Stainless steel top with stainless steel drop shelves on both sides
- The large top workspace (52.5" in length when shelves are fully extended) is ideal for holding a baby scale and other large equipment
- Key lock with two keys
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish on body, inside and out
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- 27" of vertical drawer space with fully customizable drawer configurations
- Pre-threaded mounting holes for accessories
- Industry best 5" quiet rolling casters - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Latex free cart is delivered fully assembled
- Drawer Configuration: One 6" drawer (15.2 cm) and one 21" double door storage compartment (53.3 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Body and drawers shown in Wood Vinyl Driftwood Gray.



Engineered & Manufactured in the U.S.



Three standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.



Medium Height, Medium Width Delivery Cart with Stainless Steel Top and Side Flyout Shelves, Key Lock

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300