


650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 \* Ph: (719)637-0300 Fax:(719)597-8273 \* www.harloff.com

# Sixteen Scope Storage Cabinet, with Tambour Locking Door, Drying Cabinet Package #SC16DP

The SC16DP is a state of the art stationary cabinet that stores and dries up to 16 scopes. This cabinet features two fans and a HEPA filter for drying scopes. The SC16 can be purchased with or without the HEPA filtration and drying fans. Accommodates colonoscopies, EUS Scopes, TEE probes, bronchoscopes and cystoscopes



**Specifications** 

**Construction:** 

• Sturdy, 18 gauge steel

#### Features:

- Drying fans 100-240VAC Universal Input
- HEPA filter captures up to 99.97% of particles, 0.3 microns,
- exceeding AORN scope storage recommendations
- Holds up to 16 scopes in a secure stationary location
- Back wall with padded vinyl prevents damage to scope optics
- Lockable tambour door helps protect and secure scopes
- Engineered to protect scopes and provide easy access
- Removable drip tray



Drying fan has multiple international input options

### Also Available:

Cabinet without drying package #SC16

## Color Shown is Sand Custom Colors Available

### **Assembled Cabinet Dimensions:**

- 93.25"H x 36.10"W x 18.75"D (236.9 cm H x 91.7 cm W x 47.6 cm D)
- Cabinet net weight 236 pounds (107 kg)