




650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Sixteen Scope Storage Cabinet, with Key Locking Tempered Glass Double Doors #SC16DR

The SC16DR is a state of the art stationary cabinet that stores up to 16 scopes. Accommodates colonoscopes, EUS scopes, TEE probes, bronchoscopes and cystoscopes.



**Color Show is Sand**  
**Additional Colors Available**

### Specifications

#### Construction:

- Sturdy, 18 gauge steel

#### Features:

- Holds up to 16 scopes in a secure stationary location
- Key locking tempered glass double doors
- Engineered to protect scopes and provide easy access
- Removable drip tray
- Twelve year warranty against manufacturers defect
- Adjustable hang length:
  - o High point: 82"
  - o Low point: 74"

#### Drying Package Available: #SC16DRDP Includes

- Drying fans 100-240VAC Universal Input (220 with International plug adapter # SCINTLPLUGADPTR)
- HEPA filter captures up to 99.97% of particles, 0.3 microns, exceeding AORN scope storage recommendations

#### Assembled Cart Dimensions:

- 93.25"H x 36.16"W x 20.36 "D  
(236.9 cm H x 91.8 cm W x 51.7 cm D)
- Cart net weight 304 pounds (140 kg)