

FOR ENT SCOPES

and other small scopes



SureDry
SCOPE CABINETS

**6 ENT SCOPE CAPACITY WALL MOUNT
STORAGE CABINET, KEY LOCK**

SCW2430DR

This new wall mount ENT scope cabinet allows for ENT scopes to be stored in the same safe, efficient manner of scopes in other parts of the facility

- Customizable scope holders to fit most common scope types
- Stores up to 6 ENT Scopes
- Glass paneled doors for easy viewing of scope inventory
- Air flow vents on the side of the cabinet
- Industry best 12 year limited manufacturer warranty
- Ideal for non-connected or removable video/light processor cord
- Cabinet dimensions: 30"H x 24"W x 12"D

Comes in Stainless Steel or heat sealed powder coat low carbon steel, available in over 20 color options



Engineered & Manufactured in the U.S.

This optional LED-UV disinfection strip helps to keep scopes sterilized!



UV

HARLOFF
INNOVATING STORAGE FOR CARE


UV-C disinfection range is from 185 to 280nm. Our LED runs at 270nm, close to the optimal disinfection wavelength.

View similar products online at

<https://www.harloff.com/scope-drying-and-storage-cabinets/suredry-scope-cabinets/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**

Scopes
not included



10 Reasons to Choose **HARLOFF** for Your Scope Cabinets

A Smart Choice to Meet & Exceed Current SGNA, AORN and AAMI Guidelines



1

PureDry
SCOPE CABINETS
Industry Best Positive Pressure HEPA Filtered Drying System. Designed to add Dri-Scope Aid® now or later with easy “Plug & Play” installation.

Proudly Made in the USA

6



Harloff products are proudly **Designed, Engineered and Manufactured in the United States** at our Colorado Springs, CO facility.

Optional LED UV-Disinfecting Lights

2

Help to keep your scopes sterilized with an optional UV radiation light strip as a way to disinfect your surfaces. Lights trigger in tandem with the fans after the cabinet doors are closed for scope drying.

Custom Requests - No Problem!

7

Our in-house engineering team gives us an unparalleled ability to do custom and unique projects. *Let us design the perfect scope storage system for your facility!*

Optional Wall or Cart Mounting

3

Our cabinets fit your floor plan! Use our wall-mounted cabinet for permanent installations or choose a mobile cart-mounted cabinet to move scopes from your reprocessing area directly to your treatment room.

Keyed or Keyless Locks

8

Options include latch locks, key locks or standard electronic push button locks. View

your locking options online at www.harloff.com/locking-options/



Versatile Scope Holders

4

Holds major brands including **Optim, Fuji, Olympus, Pentax, Karl Storz and others.**



Scope Cabinet Finish Options

9

Choose Stainless Steel, Heat Sealed Powder Coat in High Gloss or Textured Low Gloss, or upgrade to an **antimicrobial finish.**

12 Year Warranty

5

Our Manufacturer's Warranty is Unmatched in the industry. We build our cabinets to last and we stand behind them.



Superior Color Choices

Over 20 Colors to choose from. Cabinets have a heat sealed powder coat finish inside and out. View color options online at www.harloff.com/finish-options/

10