

FOR BRONCHOSCOPES

ENT scopes & other small scopes



ureDry
SCOPE CABINETS

**FOUR BRONCHOSCOPE OR EIGHT
ENT SCOPE CAPACITY WALL MOUNT
STORAGE CABINET, KEY LOCK**

Scopes not included

SCW2448DR


*This new wall mounted
bronchoscope cabinet
provides the same safe,
vertical drying capabilities
as that of scopes in other
parts of your facility*

- Customizable scope holders to fit most common scope types
- Holds up to 4 Bronchoscopes or up to 8 ENT Scopes
- Glass paneled doors for easy viewing of scope inventory
- Air flow vents on the side of the cabinet
- Industry best 12 year limited manufacturer warranty
- Cabinet dimensions:
48" H x 24" W x 12.75" D



**Engineered &
Manufactured in the U.S.**

**This optional LED-UV disinfection strip
helps to keep scopes sterilized!**



UV

*Comes in Stainless Steel
or heat sealed powder
coat steel, available in
over 20 color options*


HARLOFF
INNOVATING STORAGE FOR CARE

UV-C disinfection range is from
185 to 280nm. Our LED runs
at 270nm, close to the optimal
disinfection wavelength.

View similar products online at

<https://www.harloff.com/product-category/ent-and-bronchoscope-cabinets/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**



10 Reasons to Choose **HARLOFF** for Your Scope Cabinets

A Smart Choice to Meet & Exceed Current SGNA, AORN and AAMI Guidelines



1

Industry Best Positive Pressure HEPA Filtered Drying System. Designed to add Dri-Scope Aid® now or later with easy “Plug & Play” installation.

Proudly Made in the USA



Harloff products are proudly **Designed, Engineered and Manufactured in the United States** at our Colorado Springs, CO facility.

6

Optional LED UV-Disinfecting Lights

Help to keep your scopes sterilized with an optional UV radiation light strip as a way to disinfect your surfaces. Lights trigger in tandem with the fans after the cabinet doors are closed for scope drying.

Custom Requests - No Problem!

Our in-house engineering team gives us an unparalleled ability to do custom and unique projects. **Let us design the perfect scope storage system for your facility!**

Optional Wall or Cart Mounting

Our cabinets fit your floor plan! Use our wall-mounted cabinet for permanent installations or choose a mobile cart-mounted cabinet to move scopes from your reprocessing area directly to your treatment room.



Keyed or Keyless Locks

Options include latch locks, key locks or standard electronic push button locks. View your locking options online at www.harloff.com/locking-options/

Versatile Scope Holders

Holds major brands including **Optim, Fuji, Olympus, Pentax, Karl Storz and others.**

Scope Cabinet Finish Options

Choose Stainless Steel, Heat Sealed Powder Coat in High Gloss or Textured Low Gloss, or upgrade to an **antimicrobial finish.**

12 Year Warranty

Our Manufacturer's Warranty is Unmatched in the industry. We build our cabinets to last and we stand behind them.



Superior Color Choices

Over 20 Colors to choose from. Cabinets have a heat sealed powder coat finish inside and out. View color options online at www.harloff.com/finish-options/