

Standard Line Med Bin Cart with 28 – 5 inch Bins, CompX® Electronic Lock, Key Lock Override **#SL28BIN5-EKC**

Harloff medication carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.



Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Suggested Cart Option:

- 41052 - LOCKVIEW software and USB dongle with RJ-11 cable.
- Allows user to interact with electronic lock through a laptop or PC.
- Software enables mass programming, view of audit trail (who accessed the cart and when), changes in auto re-lock timing, etc.

Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 28 - 5" wide medication bins
- Dividers and labels included
- Single width narcotics drawer
- CompX® electronic lock on cabinet with BEST® brand key lock override
- Illinois lock on narcotics drawer
- Access right may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set of programmed time)
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish available in 17 colors

Left Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 7 ea – 3.25" med bin drawers with 4 – 5" bins per drawer
- 1 ea – 6.5" multipurpose drawer with adjustable dividers

Right Column Drawer Configuration:

- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 228 pounds (103.6 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 68.6 cm D)
- Approx. shipping weight: 270 pounds (122.7 kg)