

## Standard Line Med Bin Cart with 54 – 3.5 inch Bins, Key Locking #SL54BIN3

Harloff medication carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care..



### Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Construction:

- Painted steel cabinet and aluminum drawers

### Features:

- Dual column medication cart with 54 - 3.5" wide medication bins
- Dividers and labels included
- Single width narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Twelve year manufacturer's warranty

### Left Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 9 ea – 3.25" med bin drawers with 6 – 3.5" bins per drawer

### Right Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

### Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D  
(116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 235 pounds (106.8 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D  
(144.8 cm H x 119.4 cm W x 68.6 cm D)
- Approx. shipping weight: 275 pounds (125 kg)