



# MEDICATION CARTS

## SL602PC

### PUNCH CARD MEDICATION CART 750 Punch Card Capacity with a Key Lock

#### Features:

- Dual column medication cart for punch card and liquid storage with a maximum capacity of 750 punch cards
- Includes one double wide locking narcotics drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Protective gray corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers

Color shown is Cream. More than 20 additional colors are available.



#### Standard 750 Punch Card Capacity Medication Cart with BEST® Lock on Cabinet, Standard Package

Harloff medication carts provide our customers with the durability and security needed for medication storage. Standard Line med carts include narcotics storage and dividers to match your medication packaging to ensure that you receive a cart that is ready to use and will meet the demands of Long Term Care.

##### Left Column Drawer Configuration:

- 1 ea - 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawers

##### Right Column Drawer Configuration

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 9.75" locking narcotics drawer
- 2 ea - 9.75" multipurpose storage drawer

##### Assembled Cart Dimensions:

- Cart weight: 195 lbs (88.6 kg)
- 45.75" H x 43" W x 23.5" D (116.2 cm x 109.2 cm x 59.7 cm)
- Approx shipping weight: 235 pounds (106.8 kg)
- Approx Shipping Dims: 57" H x 47" W x 24" D (144.8 cm x 119.4 cm x 68.6 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Engineered & Manufactured in the U.S.**



View similar products online at  
<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300