



STC-AL6 | STC-AL8 | STC-AL10

SCOPE TRAY SYSTEM with Six, Eight or Ten Tray Capacities

Features:

- Secure scope tray holding system, with individual trays held in place by zero-slide tray retainers
- Holds multiple brands of scope transport trays
- Available workspace on the cart top for additional supplies or equipment
- 3" (7.6 cm) single wheel stem caster set, full swivel – two with brake
- Durable scratch resistant hard baked powder coat finish
- Lightweight and simple to clean
- Cart is delivered fully assembled
- Latex free
- **Scope trays not included**

Construction:

- Sturdy, all aluminum construction

Colors shown are Red and Hammertone Green. More than 20 additional colors are available.

Designed for the storage and safe transport of scope trays. Carts can be color coded to signify clean, patient ready scopes versus soiled scopes needing reprocessing. Each cart securely holds most popular brands of scope trays.



Engineered & Manufactured
in the U.S.

Optional Feature

Add a locking bar (your choice of color) to secure all trays with a single breakaway seal which can be color coded to indicate clean or soiled



STC-AL6 Dimensions:

- 37.25" H x 20.5" W x 24.375" D (94.6 cm H x 52.1 cm W x 61.9 cm D)
- Cart weight: 21.5 lbs (9.8 kg)
- Approx shipping dims: 49" H x 24" W x 29" D (124.5 cm H x 61 cm W x 73.7 cm D)
- Approx shipping weight: 52 lbs (23.6 kg)

STC-AL8 Dimensions:

- 47.25" H x 20.5" W x 24.375" D (120 cm H x 52.1 cm W x 61.9 cm D)
- Cart weight: 25.75 lbs (11.7 kg)
- Approx shipping dims: 56" H x 24" W x 29" D (142.2 cm H x 61 cm W x 73.7 cm D)
- Approx shipping weight: 61 lbs (27.7 kg)

STC-AL10 Dimensions:

- 57.25" H x 20.5" W x 24.375" D (145.4 cm H x 52.1 cm W x 61.9 cm D)
- Cart weight: 29.75 lbs (13.5 kg)
- Approx shipping dims: 49" H x 24" W x 29" D (157.5 cm H x 61 cm W x 73.7 cm D)
- Approx shipping weight: 65 lbs (29.5 kg)



View similar products online at

<https://www.harloff.com/scope-drying-and-storage-cabinets/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**