



*Breakaway bar also available
on 8 tray & 10 tray models*

STC-AL6-B

SCOPE TRANSIT SYSTEM

with Six Tray Capacity and a Breakaway Lock Bar

Features:

- Secure scope tray holding system, with individual trays held in place by zero-slide tray retainers
- Holds multiple brands of scope transport trays
- Heat sealed powder coat steel break-away lock bar (with a white antimicrobial option)
- Pack of 100 break-away seals
- Available workspace on the cart top for additional supplies or equipment
- 3" (7.6 cm) single wheel stem caster set, full swivel – two with brake
- Durable scratch resistant heat sealed powder coat finish
- Lightweight and simple to clean
- Cart is delivered fully assembled
- Latex free
- **Scope trays not included**

Construction:

- Sturdy, all aluminum construction

Colors shown are Brushed Silver and Hammertone Green. More than 20 additional colors are available.



Engineered &
Manufactured
in the U.S.

Designed for the storage and safe transport of scope trays. Carts can be color coded to signify clean, patient ready scopes versus soiled scopes needing reprocessing. Each cart securely holds most popular brands of scope trays.

STC-AL6-B Dimensions:

- 37.25" H x 20.5" W x 24.375" D
(94.6 cm H x 52.1 cm W x 61.9 cm D)
- Cart weight: 21.5 lbs (9.8 kg)
- Approx shipping dims: 49" H x 24" W x 29" D (124.5 cm H x 61 cm W x 73.7 cm D)
- Approx shipping weight: 52 lbs (23.6 kg)



Optional Feature

Breakaway seals can be color coded to indicate clean or soiled or matched to the color of your locking bar, securing all trays with a single disposable color coded lock.

View similar products online at
<https://www.harloff.com/scope-drying-and-storage-cabinets/>

HARLOFF
INNOVATING STORAGE FOR CARE

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300