



Dri-Scope Aid® STC-AL6-JET

JET CART SCOPE DRYING SYSTEM

Six Tray Capacity w/ Jet~Stream by Dri-Scope Aid®

Cart Features:

- Sturdy, all aluminum construction
- Secure scope tray holding system, with individual trays held in place by zero-slide tray retainers
- Holds up to six scope transport trays from multiple brands
- Scope hangers hold up to four scopes at once
- Eight side pockets hold disposable connectors or other supplies
- Available workspace on the cart top with a fold out shelf for extra work area, capable of holding up to two trays when extended
- 3" (7.6 cm) single wheel stem caster set, full swivel – two with brake
- Durable scratch resistant hard baked powder coat finish
- Lightweight and simple to clean
- Cart is delivered fully assembled
- Latex free
- **Scopes and scope trays not included**

Jet~Stream Features:

- Includes the Jet~Stream connector storage system from Dri-Scope Aid®
- Provides compact drying for up to two scopes at once
- Significantly improves the drying of internal channels
- HEPA filter offers a 99.97% retention rate of particles that have a size > 0.3µm
- Provides constant air delivery through all internal endoscope channels
- Short drying cycle fulfills forced air recommendations and provides a quick turn-around time
- Self-adjusting high-pressure regulation
- Hands-free operation with visual display of purging times
- Simple to attach individual scope connectors
- Works with all major endoscope brands
- Disposable connectors meet SGNA and AAMI guidelines



Color shown is White.
More than 20 additional
colors are available.

Jet~Stream device
sold separately

Dimensions:

- 37.25" H x 20.5" W x 24.375" D cart only (94.6 cm H x 52.1 cm W x 61.9 cm D)
- 74" H x 43" W x 31" D fully extended (188 cm H x 109.2 cm W x 78.7 cm D)
- Cart weight: 51 lbs (23.1 kg)
- Approx shipping dims: 49" H x 24" W x 29" D (124.5 cm H x 61 cm W x 73.7 cm D)
- Approx shipping weight: 81 lbs (36.7 kg)



Engineered &
Manufactured
in the U.S.

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/scope-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300