


# STC-SH4L | STC-SH6L STC-SH8L | STC-SH10L

## SCOPE TRANSIT SYSTEM

with Four, Six, Eight or Ten Tray Capacities



### Features:

- Sturdy, steel construction
- Secure scope tray holding system, with individual trays held in place by zero-slide tray retainers
- Holds multiple brands of scope transport trays, including standard (20.625" W x 16.75" D x 4.125" H) trays and larger Advantage Plus™ and Cleanascope™ style trays (20.375" W x 18.25" D x 4.75" H)
- Available workspace on the cart top for additional supplies or equipment
- 3" casters – two with brake
- Durable scratch resistant hard baked powder coat or stainless steel finish
- Lightweight and simple to clean
- Latex free cart is delivered fully assembled
- **Scope trays not included**

Designed for the storage and safe transport of scope trays. Carts can be color coded to signify clean, patient ready scopes versus soiled scopes needing reprocessing.



Engineered & Manufactured  
in the U.S.

Colors shown are Beige, Hammer Black, Red and HT Green. More than 20 additional colors or stainless steel finishes are available.

STC-SH10L

STC-SH6L

STC-SH8L

STC-SH4L

#### STC-SH4L Dimensions:

- 28.63" H x 24.25" W x 23.88" D  
(72.7 cm H x 61.6 cm W x 60.7 cm D)
- Cart weight: 60.3 lbs (27.4kg)

#### STC-SH6L Dimensions:

- 39.38" H x 24.25" W x 24.63" D  
(100 cm H x 61.6 cm W x 62.6 cm D)
- Cart weight: 75 lbs (34 kg)

#### STC-SH8L Dimensions:

- 50.38" H x 24.13" W x 24.63" D  
(128 cm H x 61.3 cm W x 62.6 cm D)
- Cart weight: 90 lbs (40.8 kg)

#### STC-SH10L Dimensions:

- 61.38" H x 24.13" W x 24.63" D  
(155.9 cm H x 61.3 cm W x 62.6 cm D)
- Cart weight: 105 lbs (47.6 kg)

Add an optional locking bar (your choice of color) to secure all trays with a single breakaway seal which can be color coded to indicate clean or soiled



View similar products online at

<https://www.harloff.com/scope-drying-and-storage-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300