



STC-AL6 | STC-AL8 | STC-AL10

SCOPE TRANSIT SYSTEM OPTIONS

Breakaway locks w/color coded seals or cart bodies

Change the appearance, expand the functionality and broaden the usefulness of your Harloff cart by purchasing optional built-in add-on products.

Breakaway Lock Bar:

- Secure your scope trays with a breakaway lock bar
- Comes with a pack of 100 breakaway seals
- Enhances identification of trays allowing you to more easily keep track of which scope trays are clean or soiled
- Each breakaway seal has a unique identifying number, allowing you to track which staff is responsible for tray reprocessing
- Helps prevent tampering by restricting tray access
- Comes in a heat sealed powder coat steel finish in over twenty color choices including an antimicrobial option

Color Coded Seals:

- Breakaway seals come in a default yellow, but can be purchased in several color options
- Color code your seals: green for clean, red for soiled or any other color options to match your procedures
- Adds an additional identifying feature



Color Coded Cart Bodies:

- Color code your cart body: green for clean, red for soiled or any other color options to match your procedures
- Color code to specific reprocessing areas
- Our durable heat sealed powder coat aluminum finish and open design is easy to clean and sterilize



Engineered & Manufactured in the U.S.

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/scope-drying-and-storage-cabinets/>

Toll Free: (800) 433-4064
medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300

