



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## V-Series, Three Drawer Procedure Cart with Key Lock #V24-3K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

### Specifications

#### Construction:

- Steel construction

#### Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 24" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

#### Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

#### Assembled Cart Dimensions

- 34" H x 29.5" W x 22" D  
(86.4 cm H x 74.9 cm W x 55.9 cm D)
- Cart weight: 88 pounds (39.9 kg)
- Approximate shipping weight: 110 pounds (49.9 kg)

#### Drawer Configuration:

- (1) – 6" (15.2 cm)
  - (2) – 9" (22.9 cm)
- Integrated Pull Handle on all drawers



### Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>