



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

V-Series, Two Drawer Procedure Cart with Locking Storage Compartment #V30-2K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

Specifications



Construction:

- Steel construction

Features:

- Key locking drawers and storage compartment; includes four keys
- Aluminum side rail / push handle
- Tall cabinet
- 21" storage compartment with locking doors and adjustable shelf
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6cm H x 74.9 W x 55.9 D)
- Cart weight: 90 lbs (40.8 kg)
- Approx shipping weight: 113 pounds (51.2 kg)

Drawer Configuration:

- (1) – 3" (7.6 cm)
 - (1) – 6" (15.2 cm)
 - Storage compartment – 21" (53.3 cm)
- Integrated Pull Handle on all drawers

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>