



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## V-Series, Six Drawer Crash Cart, Deluxe Emergency Accessory Package #V30-6EMG2

This high capacity crash cart is designed for healthcare facilities under budget pressure requiring a durable, fully equipped emergency cart offering Harloff's quality construction at a value price.

This emergency cart includes all commonly required accessories and is perfect for hospitals, clinics, dialysis centers, doctors' offices, urgent care facilities and nursing homes needing to buy carts to store, transport and dispense critical supplies used for emergency situations.



### Specifications

#### Construction:

- Steel construction

#### Features:

- Break-away lock bar
- 100 break-away seals included
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- Minimal assembly required at delivery
- Latex free
- Durable, hard-baked powder coat finish

#### DELUXE Package Accessories Include:

- Plastic top
- Cardiac board and brackets
- Powder coated steel top rail / push handle
- Side mounted drop shelf
- IV pole
- Oxygen tank holder (D/E size cylinders)
- (3) 6-Compartment Drawer Trays
- (3) 7-Compartment Drawer Trays

#### Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

Integrated Pull Handle on all drawers

#### Sixteen available color choices:

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

#### Assembled Cart Dimensions

- 67.4" H x 33.25" W x 22" D (171.2 cm H x 84.5 W x 55.9 D) (height includes IV Pole)
- Cart weight: 138 lbs (62.6 kg)
- Approx shipping weight: 157 pounds (71.2 kg)