



MEDICATION CARTS

VLT240PC

VALUE PUNCH CARD MED CART 300 Punch Card Capacity with a Key Lock

Features:

- Single column medication cart for punch card and liquid storage with a maximum capacity of 300 punch cards
- Includes one in-drawer single width locking narcotics box in top 9" drawer; comes with an Illinois lock with four keys
- Tall cabinet with key locking system with four keys
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Taller drawer height for increased capacity
- Ball-bearing, full extension drawer slides
- Painted steel top rail/push handle
- Antimicrobial plastic top, removable
- Protective gray corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Adjustable punch card row dividers
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Certified "Latex Free"
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered &
Manufactured
in the U.S.

Color shown is Cream. More than
20 additional colors are available.



Value 300 Punch Card Capacity Medication Cart with Standard Key Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 1 - ea 9.75" drawer with 1 divider and locking narcotics box
- 2 - ea 9.75" drawers with two adjustable dividers

Assembled Cart Dimensions:

- Writing surface: 43.25" (109.9 cm)
- 47" H x 21.3" W x 21.5" D (119.4 cm H x 54.1 cm W x 54.6 cm D)
- Cart weight: 110 pounds (50 kg)
- Approx shipping weight: 155 pounds (70.45 kg)
- Approx shipping dims: 48"H x 33" W x 26" D (121.9 cm x 83.8 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at

<https://www.harloff.com/product-category/punch-card-carts-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300