

Value Line Med Bin Cart with 20 – 5 inch Bins, Basic Electronic Pushbutton Lock with Key Lock Override #VLT24EBIN5

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with basic electronic lock system, key lock override with 4 keys
- Factory set user and supervisor code
- Lock can be programmed in field
- Internal narcotics box in bottom 9" drawer
- Pull out side shelf
- Includes 24 Harloff 5" medication bins with dividers and labels
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Antimicrobial plastic top, removable
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawers
- 6 - ea 3" bin drawers including 4 - 5" bins each
- 1 - ea 9" drawer with 1 divider and locking narcotics box

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 165 pounds (75 kg)
- Approx. shipping weight: 210 pounds (95.5kg)
- Approx. shipping dims: 57"H x 33" W x 26" D (144.8cm H x 83.8cm W x 66cm D)

Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand