

## Value Line Punch Card Medication Cart with Basic E-Lock and Key Lock Override #VLT480EPC

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



### Construction:

- Painted steel cabinet and drawers

### Features:

- Tall cabinet with basic electronic pushbutton lock, key lock override with 4 keys
- Factory set user and supervisor code
- Lock can be programmed in the field
- Capacity of up to 480 punch cards
- Internal locking narcotics box in top 9" drawer
- Painted steel pull-out side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel, - one with brake, one directional
- Painted steel top rail/push handle
- Antimicrobial plastic top, removable
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

### Internal Drawer Dimensions:

- 27.5" W x 17" D (68.75 cm W x 43.2 cm D)

### Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 1 - ea 9.75" drawer with 1 divider and locking narcotics box
- 2 - ea 9.75" drawers with two adjustable dividers

### Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D  
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 155 pounds (70.45 kg)
- Approx. shipping weight: 190 pounds (89 kg)
- Approx. shipping dims: 57" H x 37" W x 27" D  
(144.8cm H x 93.9cm W x 68.6cm D)

### Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand