



MEDICATION CARTS

VLT480PC

PUNCH CARD MEDICATION CART 600 Punch Card Capacity with a Key Lock

Features:

- Single column medication cart for punch card and liquid storage with a maximum capacity of 600 punch cards
- Includes one single width locking narcotics drawer with an Illinois lock
- Key lock on cabinet with four keys
- Painted steel cabinet with aluminum drawers
- Durable scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Painted steel top rail / push handle
- Removeable plastic top
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Protective gray corner bumpers
- Adjustable punch card row dividers
- Body has pre-threaded accessory mounting holes
- Cart is delivered fully assembled
- Latex free



Color shown is Cream. More than 20 additional colors are available.

Value 600 Punch Card Capacity Medication Cart with Key Lock on Cabinet, Standard Package

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget. A complete line of accessories and narcotics storage optimize your cart's usefulness to nursing staff.

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 1 - ea 9.75" drawer with one divider and locking narcotics box
- 2 - ea 9.75" drawers with two adjustable dividers

Assembled Cart Dimensions:

- Cart weight: 155 pounds (70.45 kg)
- 46.75" H x 33" W x 22.5" D (121.3 cm H x 82.5 cm W x 57.2 cm D)
- Approx. shipping weight: 185 pounds (89 kg)
- Approx Shipping dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.



View similar products online at <https://www.harloff.com/product-category/punch-card-carts-cabinets/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**